


Figure S7. Recommended protocol and technical aspects for non-adherent cells

NON-ADHERENT CELLS



Metabolome coverage

ESI	HILIC (Log P < 1)	RP modified (-3 < Log P < 7)	RP (Log P > -1)
+	379 metabolites: heterocyclic comp., aas, carbohydrates, benzenoids, acylcarnitines...	1089 metabolites: heterocyclic comp., benzenoids, acylcarnitines, steroids, phospholipids, aas, amines...	1001 metabolites: heterocyclic comp., benzenoids, steroids, aas, organic oxygen...
-	110 metabolites: carbohydrates, FAs, organic acids...	268 metabolites: FAs, ecosanoids, phospholipids, carbohydrates, benzenoids, organic acids, flavonoids...	241 metabolites: FAs, ecosanoids, phospholipids, benzenoids, flavonoids, organic acids...
+/-	486 metabolites: aas, heterocyclic comp., nucleotides, carbohydrates, benzenoids...	973 metabolites: flavonoids, heterocyclic comp., FAs, sterpoids...	818 metabolites: FAs, steroids, BAs, flavonoids, benzenoids, aaa...

aas: amino acids; FAs: fatty acids; BAs: bile acids

Tips

- Ideally adjust extraction volume to cell number or packed cell volume, plan ahead to have 1 replicate per condition to this end.
- Perform the optional wash step only if contamination from media metabolites is relevant.
- Once the media is removed work as quickly as possible to prevent alterations in the metabolome, ideally process a max of 6 samples at a time until the step where the extraction solvent is added.
- Use ice-cold extraction solvent, mix well and keep the samples on ice to minimize enzymatic metabolite degradation.
- Centrifuge at least 2x at max speed at 4 °C to ensure that no cell debris is transferred to the HPLC vial.

[*] Recommended conditions based on our experience